

bsm-63210R

[Primary Antibody]

TTR/Prealbumin Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 5G6

GeneID: 7276

Target: TTR/Prealbumin

Immunogen: A synthesized peptide derived from human Transthyretin: 100-147.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Applications: **WB** (1:1000-2000)

IHC-P (1:100-200)

IHC-F (1:100-200)

IF (1:50-200)

Flow-Cyt (1:50-100)

ICC/IF (1:50-200)

IP (1:20-50)

Reactivity: Human

Predicted MW.: 16 kDa

Subcellular Location: Secreted ,Cytoplasm