
NFKB2 Recombinant Rabbit mAb

Catalog Number: bsm-63222R

Target Protein: NFKB2

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: 5A11

Isotype: IgG

Applications: WB (1:500-2000), ICC/IF (1:50-200), IP (1:20-50)

Reactivity: Human, Mouse, Rat

Predicted MW: 97 kDa

Detected MW: 120 kDa

Entrez Gene: 4791

Swiss Prot: Q00653

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: NFkB-p100 a transcription factor of the nuclear factor-kappaB (NFkB) group. Precursor of the p52 subunit of the nuclear factor NF-kappa-B, which binds to the kappa-B consensus sequence 5'-GGRNNYCC-3', located in the enhancer region of genes involved in immune response and acute phase reactions.