

**bsm-63222R****[ Primary Antibody ]****NFKB2 Recombinant Rabbit mAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**— DATASHEET —****Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**CloneNo.:** 5A11**GeneID:** 4791**SWISS:** Q00653**Target:** NFKB2**Immunogen:** A synthesized peptide derived from human NF- $\kappa$ B p52: 700-900/900.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** NF $\kappa$ B-p100 a transcription factor of the nuclear factor-kappaB (NF $\kappa$ B) group. Precursor of the p52 subunit of the nuclear factor NF-kappa-B, which binds to the kappa-B consensus sequence 5'-GGRNNYYCC-3', located in the enhancer region of genes involved in immune response and acute phase reactions.**Applications:** **WB** (1:500-2000)**ICC/IF** (1:50-200)**IP** (1:20-50)**Reactivity:** Human, Mouse, Rat**Predicted  
MW.:** 97 kDa**Subcellular  
Location:** Nucleus