

bsm-63253R

[Primary Antibody]

MYB Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 5H4

GeneID: 4602

SWISS: P10242

Target: MYB

Immunogen: A synthesized peptide derived from human c Myb: 600-640/640.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

Applications: **WB** (1:1000-2000)

IHC-P (1:100-200)

IHC-F (1:100-200)

IF (1:100-200)

Reactivity: Human

**Predicted
MW.:** 72 kDa

**Subcellular
Location:** Nucleus