

**bsm-63253R**

**[ Primary Antibody ]**

## MYB Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**CloneNo.:** 5H4

**GeneID:** 4602

**SWISS:** P10242

**Target:** MYB

**Immunogen:** A synthesized peptide derived from human c Myb: 600-640/640.

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

**Applications:** **WB** (1:1000-1:2000)

**IHC-P** (1:100-1:200)

**IHC-F** (1:100-1:200)

**IF** (1:100-1:200)

**Reactivity:** Human

**Predicted  
MW.:** 72 kDa

**Subcellular  
Location:** Nucleus