

**bsm-63258R****[ Primary Antibody ]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**EGFR(L858R mutated) Recombinant Rabbit mAb****— DATASHEET —**

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>WB</b> (1:1000-1:2000) <b>IHC-P</b> (1:100-1:200) <b>IHC-F</b> (1:100-1:200) <b>IF</b> (1:50-1:200) <b>ICC/IF</b> (1:50-1:200)  <b>Reactivity:</b> Human   <b>Predicted MW.:</b> 134 kDa  <b>Subcellular Location:</b> Secreted ,Cell membrane
<b>Clonality:</b> Recombinant	<b>CloneNo.:</b> 5C5	
<b>GeneID:</b> 1956	<b>SWISS:</b> P00533	
<b>Target:</b> EGFR(L858R mutated)		
<b>Immunogen:</b> A synthesized peptide derived from human EGFR (L858R mutated): FG-R-AK.		
<b>Purification:</b> affinity purified by Protein A		
<b>Storage:</b> 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> Receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. Known ligands include EGF, TGFA/TGF-alpha, AREG, epigen/EPGN, BTC/betacellulin, epiregulin/EREG and HBEGF/heparin-binding EGF.		