

bsm-63286R

[Primary Antibody]

IL-15 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) Reactivity: Human Predicted MW.: 18 kDa Subcellular Location: Secreted ,Cell membrane ,Cytoplasm ,Nucleus
Clonality: Recombinant	CloneNo.: 5G3	
GeneID: 3600	SWISS: P40933	
Target: IL-15		
Immunogen: A synthesized peptide derived from human IL 15: 100-162.		
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: Cytokine that stimulates the proliferation of T-lymphocytes. Stimulation by IL-15 requires interaction of IL-15 with components of IL-2R, including IL-2R beta and probably IL-2R gamma but not IL-2R alpha.		