

bsm-63290R

[Primary Antibody]

phospho-ERBB3 (Y1289) Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-1:2000) Reactivity: Human Predicted MW.: 148 kDa Subcellular Location: Secreted ,Cell membrane
Clonality: Recombinant	CloneNo.: 5D3	
GeneID: 2065	SWISS: P21860	
Target: phospho-ERBB3 (Y1289)		
Immunogen: A synthesized peptide derived from human ErbB3 around the phosphorylation site of Y1289: EQG-pY-EEM.		
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: Binds and is activated by neuregulins and NTAK.		