

**bsm-63291R**

**[ Primary Antibody ]**

## **BMP15 Recombinant Rabbit mAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **— DATASHEET —**

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**CloneNo.:** 5D4

**GeneID:** 9210

**SWISS:** O95972

**Target:** BMP15

**Immunogen:** A synthesized peptide derived from human BMP15: 350-392/392.

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** May be involved in follicular development. Oocyte-specific growth/differentiation factor that stimulates folliculogenesis and granulosa cell (GC) growth.

**Applications:** **WB** (1:1000-1:2000)

**Flow-Cyt** (1:20-1:100)

**Reactivity:** Human, Mouse, Rat

**Predicted  
MW.:** 45 kDa

**Subcellular  
Location:** Secreted