

bsm-63293R

[Primary Antibody]

phospho-ERBB3 (Y1222) Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 6E2

GeneID: 2065

SWISS: P21860

Target: ERBB3 (Y1222)

Immunogen: A synthesized peptide derived from human ErbB3 around the phosphorylation site of Y1222: LEELG-pY-EY.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Binds and is activated by neuregulins and NTAK.

Applications: **WB** (1:1000-2000)

ICC/IF (1:50-200)

IP (1:20-50)

Reactivity: Human

Predicted MW.: 148 kDa

Subcellular Location: Secreted ,Cell membrane