

bsm-63294R

[Primary Antibody]

phospho-Alpha synuclein (Y125) Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-1:2000) Reactivity: Human Predicted MW.: 14 kDa Subcellular Location: Cell membrane ,Cytoplasm ,Nucleus
Clonality: Recombinant	CloneNo.: 2E7	
GeneID: 6622	SWISS: P37840	
Target: phospho-Alpha synuclein (Y125)		
Immunogen: A synthesized peptide derived from human alpha Synuclein around the phosphorylation site of Y125: NEA-pY-EM.		
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: Neuronal protein that plays several roles in synaptic activity such as regulation of synaptic vesicle trafficking and subsequent neurotransmitter release.		