

bsm-63297R

[Primary Antibody]

IL-9 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 1B5

GeneID: 3578

SWISS: P15248

Target: IL-9

Immunogen: A synthesized peptide derived from human IL 9: 100-144.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Multifunctional cytokine secreted mainly by T-helper 2 lymphocytes and also mast cells or NKT cells that plays important roles in the immune response against parasites. Affects intestinal epithelial permeability and adaptive immunity.

Applications: WB (1:1000-2000)

Reactivity: Human

Predicted MW.: 16 kDa

Subcellular Location: Secreted