

**bsm-63297R****[ Primary Antibody ]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**IL-9 Recombinant Rabbit mAb****— DATASHEET —**

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> WB (1:1000-1:2000)  <b>Reactivity:</b> Human  <b>Predicted MW.:</b> 16 kDa  <b>Subcellular Location:</b> Secreted
<b>Clonality:</b> Recombinant	<b>CloneNo.:</b> 1B5	
<b>GeneID:</b> 3578	<b>SWISS:</b> P15248	
<b>Target:</b> IL-9		
<b>Immunogen:</b> A synthesized peptide derived from human IL 9: 100-144.		
<b>Purification:</b> affinity purified by Protein A		
<b>Storage:</b> 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> Multifunctional cytokine secreted mainly by T-helper 2 lymphocytes and also mast cells or NKT cells that plays important roles in the immune response against parasites. Affects intestinal epithelial permeability and adaptive immunity.		