

bsm-63303R

[Primary Antibody]

CTNNBIP1 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 12C4

GeneID: 56998

SWISS: Q9NSA3

Target: CTNNBIP1

Immunogen: A synthesized peptide derived from human beta Catenin interacting protein 1: 1-30.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Prevents the interaction between CTNNB1 and TCF family members, and acts as a negative regulator of the Wnt signaling pathway.

Applications: **WB** (1:1000-1:2000)

Flow-Cyt (1:20-1:100)

IP (1:20-1:50)

Reactivity: Human

**Predicted
MW.:** 9 kDa

**Subcellular
Location:** Cytoplasm ,Nucleus