

bsm-63305R**[Primary Antibody]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

phospho-HER2 (Y1139) Recombinant Rabbit mAb**— DATASHEET —**

| | | |
|--|----------------------|--|
| Host: Rabbit | Isotype: IgG | Applications: WB (1:1000-2000) IP (1:20-50) Reactivity: Human Predicted MW.: 138 kDa Subcellular Location: Cell membrane ,Cytoplasm ,Nucleus |
| Clonality: Recombinant | CloneNo.: 1H2 | |
| GeneID: 2064 | SWISS: P04626 | |
| Target: HER2 (Y1139) | | |
| Immunogen: A synthesized peptide derived from human ErbB2 around the phosphorylation site of Y1139: PE-pY-VN. | | |
| Purification: affinity purified by Protein A | | |
| Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles. | | |
| Background: Protein tyrosine kinase that is part of several cell surface receptor complexes, but that apparently needs a coreceptor for ligand binding. Essential component of a neuregulin-receptor complex, although neuregulins do not interact with it alone. | | |