

bsm-63305R**[Primary Antibody]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

phospho-HER2 (Y1139) Recombinant Rabbit mAb**— DATASHEET —****Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**CloneNo.:** 1H2**GeneID:** 2064**SWISS:** P04626**Target:** phospho-HER2 (Y1139)**Immunogen:** A synthesized peptide derived from human ErbB2 around the phosphorylation site of Y1139: PE-pY-VN.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** Protein tyrosine kinase that is part of several cell surface receptor complexes, but that apparently needs a coreceptor for ligand binding. Essential component of a neuregulin-receptor complex, although neuregulins do not interact with it alone.**Applications:** **WB** (1:1000-1:2000)
IP (1:20-1:50)**Reactivity:** Human**Predicted
MW.:** 138 kDa**Subcellular
Location:** Cell membrane ,Cytoplasm
Nucleus