

bsm-63307R

[Primary Antibody]

IL12B Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 3A9

Target: IL12B

Immunogen: A synthesized peptide derived from human IL 12B: 281-328.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Applications: **WB** (1:1000-1:2000)

Flow-Cyt (1:20-1:100)

ICC/IF (1:50-1:200)

Reactivity: Human, Mouse, Rat

Predicted MW.: 37 kDa