

**bsm-63309R**

**[ Primary Antibody ]**

## UBA52 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**CloneNo.:** 10A2

**GeneID:** 7311

**Target:** UBA52

**Immunogen:** A synthesized peptide derived from human UBA52: 70-97.

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** UBA52 is a fusion protein consisting of ubiquitin at the N terminus and ribosomal protein L40 at the C terminus, a C-terminal extension protein (CEP). Multiple processed pseudogenes derived from this gene are present in the genome.

**Applications:** **WB** (1:1000-2000)

**IHC-P** (1:100-200)

**IHC-F** (1:100-200)

**IF** (1:50-200)

**Flow-Cyt** (1:50-100)

**ICC/IF** (1:50-200)

**Reactivity:** Human, Mouse, Rat

**Predicted  
MW.:** 15 kDa

**Subcellular  
Location:** Cytoplasm ,Nucleus