

bsm-63309R

[Primary Antibody]

UBA52 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 10A2

GeneID: 7311

Target: UBA52

Immunogen: A synthesized peptide derived from human UBA52: 70-97.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: UBA52 is a fusion protein consisting of ubiquitin at the N terminus and ribosomal protein L40 at the C terminus, a C-terminal extension protein (CEP). Multiple processed pseudogenes derived from this gene are present in the genome.

Applications: **WB** (1:1000-1:2000)

IHC-P (1:100-1:200)

IHC-F (1:100-1:200)

IF (1:50-1:200)

Flow-Cyt (1:20-1:100)

ICC/IF (1:50-1:200)

Reactivity: Human, Mouse, Rat

**Predicted
MW.:** 15 kDa

**Subcellular
Location:** Cytoplasm ,Nucleus