

AbBy Fluor® 633 酪胺 TSA AbBy Fluor® 633

产品编号: IHC049

保存条件: -20℃避光保存。

产品介绍:

Bioss SKU	Bioss Name	激发峰波长Ex. (nm)/发射峰波长Em. (nm)
IHC041	TSA AbBy Fluor® 405	405/420
IHC042	TSA AbBy Fluor® 350	362/436
IHC043	TSA AbBy Fluor® 488	488/519
IHC044	TSA AbBy Fluor® 514	517/539
IHC045	TSA AbBy Fluor® 532	534/553
IHC046	TSA AbBy Fluor® 555	555/570
IHC047	TSA AbBy Fluor® 568	579/603
IHC048	TSA AbBy Fluor® 594	594/615
IHC049	TSA AbBy Fluor® 633	633/648
IHC050	TSA AbBy Fluor® 647	650/670
IHC051	TSA AbBy Fluor® 680	681/704
IHC052	TSA Cy7	750/780