
TDI-011536

产品编号: D51764

CAS: 2687970-96-1

分子式: C₁₉H₁₆N₄O₂S

纯度: 98.7%

InChi Key: InChIKey=LJNPGPUDPDQGLE-QOCHGBHMSA-N

Smiles: C(N1\C(=N\C(=O)C=2C=3C(NC2)=NC=CC3)\SC=C1CO)C4=CC=CC=C4

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 100 mg/mL(274.41 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TDI-011536 is a potent Lats kinase inhibitor. TDI-011536 interrupts Hippo-Yap signaling and initiates the proliferation of lesioned heartmuscle cells. TDI-011536 can be used in studies of organ conservation and regeneration.