

---

## ITMN4077

产品编号: D51769

CAS: 790305-05-4

分子式: C<sub>26</sub>H<sub>40</sub>N<sub>4</sub>O<sub>8</sub>S

纯度: 98.85%

InChi Key: InChIKey=MXGKXZHTCMUVEE-JXUUZEIMSA-N

Smiles: CC(C)(C)OC(N[C@@H](CCCCC/C=C/[C@H](C1)[C@]1(C(NS(C1CC1)(=O)=O)=O)NC([C@H](C1)N2C[C@@H]1O)=O)C2=O)=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: Antiviral

溶解性: DMSO: 55 mg/ml (96.72 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ITMN4077 has antiviral activity against Hepatitis C virus infected in human HuH7 cells assessed as inhibition of HCV subgenomic replicon replication with an EC<sub>50</sub> of 2.131μM.