

---

## PCO 400

产品编号: D51771

CAS: 121055-10-5

分子式: C<sub>17</sub>H<sub>17</sub>NO<sub>4</sub>

纯度: 98.06%

InChi Key: InChIKey=WJBHHTPFTVKZCV-CVEARBPZSA-N

Smiles: O([C@@H]1C=2C(OC(C)(C)[C@H]1O)=CC=C(C#N)C2)C3=CC(=O)CC3

外观: 固体粉末

作用通路: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PCO 400 is an analog of comakalim that acts as a selective and potent ATP-sensitive K<sup>+</sup> channel opener.