
Soquinolol

产品编号: D51774

CAS: 61563-18-6

分子式: C₁₇H₂₆N₂O₃

纯度: 98.59%

InChi Key: InChIKey=GMJHUSJLZXFFQJ-UHFFFAOYSA-N

Smiles: O=CN1CC=2C=CC=C(OCC(O)CNC(C)(C)C)C2CC1

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Adrenergic Receptor

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Soquinolol is a beta-adrenergic receptor antagonist with weak local anesthetic activity. Soquinolol has good enteral efficacy and a long duration of action.