
Wnt pathway activator 2

产品编号: D51914

CAS: 1360540-82-4

分子式: C₁₇H₁₅NO₄

纯度: 98.28%

InChi Key: InChIKey=DMBCGEJWBLNDP-UHFFFAOYSA-N

Smiles: O=C(C1=NC=CC=C1)CCC(=O)C2=CC=C3OCCOC3=C2

外观: 固体粉末

作用通路: Cytoskeletal Signaling; Stem Cells

作用靶点: Wnt/beta-catenin

溶解性: DMSO : 90.0 mg/mL (302.7 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Wnt pathway activator 2 is a potent Wnt activator with an IC₅₀ of 13 nM.