
BiP inducer X

产品编号: D51916

CAS: 101714-41-4

分子式: C₉H₇NO₃S

纯度: 98%

InChi Key: InChIKey=SVFLBLCWKKQKDW-UHFFFAOYSA-N

Smiles: N#CCCC(c(cc1)cc(O)c1O)=O

外观: 固体粉末

作用通路: Apoptosis

作用靶点: Apoptosis

溶解性: DMSO : 45 mg/mL (215.08 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BiP inducer X is a selective inducer of immunoglobulin heavy chain binding protein (BiP)/GRP78 and an ER chaperone inducer. BiP inducer X prevents cell death in neurons and retinal cell lines.