
IWP L6

产品编号: D51955

CAS: 1427782-89-5

分子式: C₂₅H₂₀N₄O₂S₂

纯度: 99.11%

InChi Key: InChIKey=QESQGTFWEQMCMH-UHFFFAOYSA-N

Smiles: O=C(CSc1nc2CCSc2c(=O)n1-c1cccc1)Nc1ccc(cn1)-c1cccc1

外观: 固体粉末

作用通路: Stem Cells

作用靶点: Porcupine

溶解性: DMSO:47.3 mg/mL (100 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: IWP-L6 is an extremely effective Porcn inhibitor (EC₅₀: 0.5 nM).