
GSK467

产品编号: D52010

CAS: 1628332-52-4

分子式: C₁₇H₁₃N₅O₂

纯度: 100%

InChi Key: InChIKey=ZTYRLXUTLYBVHH-UHFFFAOYSA-N

Smiles: O=c1[nH]c(Oc2cnn(Cc3ccccc3)c2)nc2cnccc12

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Demethylase; Histone Methyltransferase

溶解性: DMSO:20.83 mg/mL (65.23 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK467 is a potent and selective inhibitor of KDM5 (JARID1)(K_i: 10 nM).