
F-15599

产品编号: D52044

CAS: 635323-95-4

分子式: C₁₉H₂₁ClF₂N₄O

纯度: 98%

InChi Key: InChIKey=WAAXKNFGOFTGLP-UHFFFAOYSA-N

Smiles: Cc1cnc(CNCC2(F)CCN(CC2)C(=O)c2ccc(F)c(Cl)c2)nc1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 22.5 mg/mL (56.98 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: F-15599 is a highly selective G-protein biased 5-HT_{1A} receptor agonist(Ki : 3.4 nM).