
ML339

产品编号: D52087

CAS: 2579689-83-9

分子式: C₂₆H₃₂ClN₃O₅

纯度: 98%

InChi Key: InChIKey=SSPYAPRDKNCABY-LDLYASANSA-N

Smiles: O=C(NC1=CC=CC=C1Cl)CN2[C@H]3C[C@@H](NC(C4=CC(OC)=C(OC)C(OC)=C4)=O)C[C@@H]2CCC3

外观: 固体粉末

作用通路: Autophagy; GPCR/G Protein; Immunology/Inflammation

作用靶点: CXCR

溶解性: DMSO : 125 mg/mL (249.00 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ML339 a selective inhibitor of CXCR6(IC₅₀ = 140 nM) with no response when screened against CXCR5 and CXCR4.