
LY3009120

产品编号: D52108

CAS: 1454682-72-4

分子式: C₂₃H₂₉FN₆O

纯度: 99.54%

InChi Key: InChIKey=HHCBMISMPSAZBF-UHFFFAOYSA-N

Smiles: CNc1ncc2cc(c(C)nc2n1)-c1cc(NC(=O)NCCC(C)(C)C)c(F)cc1C

外观: 固体粉末

作用通路: Autophagy; MAPK

作用靶点: Autophagy; Raf

溶解性: DMSO: 10 mg/mL (23.56 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LY03009120 is a potent pan-Raf inhibitor with IC₅₀ of 44 nM, 31-47 nM, and 42 nM for A-raf, B-Raf, and C-Raf in A375 cells, respectively. Phase 1.