
BMS-935177

产品编号: D52109

CAS: 1231889-53-4

分子式: C₃₁H₂₆N₄O₃

纯度: 99.89%

InChi Key: InChIKey=TVJRDCQUZMGBAB-UHFFFAOYSA-N

Smiles: O=C(N)C1=CC=C(C=2C=CC=C(C2C)N3C=NC=4C=CC=CC4C3=O)C=5C6=CC=C(C=C6NC15)C(O)(C)C

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: BTK

溶解性: DMSO: 117mg/mL (232.8mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BMS-935177 is a reversible BTK inhibitor with IC₅₀ value of 3 nM.