
3-Cyano-7-hydroxycoumarin

产品编号: D52120

CAS: 19088-73-4

分子式: C₁₀H₅NO₃

纯度: 98.03%

InChi Key: InChIKey=IJQYTHQDUDCJEQ-UHFFFAOYSA-N

Smiles: N#CC1=CC=CC(=O)=CC2OC1=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 3-Cyano-7-hydroxycoumarin is a macrophage inhibitor of macrophage inhibitory factor (MIF) interconjugate isomerase with a K_i of 2.9 μ M. 3-Cyano-7-hydroxycoumarin belongs to the group of coumarin derivatives, which can be used as molecular probes and fluorescent dyes.