
Calcium Channel antagonist 4

产品编号: D52161

CAS: 687574-49-8

分子式: C₂₃H₂₆N₂O₄S

纯度: 98.09%

InChi Key: InChIKey=DVJRGJOPKYGOOQ-UHFFFAOYSA-N

Smiles: O=C(NCC=1C=CC=CC1)C=2OC=3C=CC(=CC3C2C)S(=O)(=O)N4CCC(C)CC4

外观: 固体粉末

作用通路: Membrane transporter/Ion channel;Metabolism

作用靶点: Calcium Channel

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Calcium Channel antagonist 4 is a voltage-gated calcium channel inhibitor (IC₅₀ : 5-20 μM).