
c-Kit-IN-3

产品编号: D52200

CAS: 2363169-01-9

分子式: C₂₆H₂₀ClF₃N₂O₄

纯度: 99.43%

InChi Key: InChIKey=YOHLRCOPRAVUCJ-UHFFFAOYSA-N

Smiles: O=C(NC1=CC=C(OC=2C=CN=C3C=C(OC)C(OC)=CC32)C=C1)CC4=CC=C(Cl)C(=C4)C(F)(F)F

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Kit

溶解性: DMSO: 100 mg/mL (193.46 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: c-Kit-IN-3 is a selective inhibitor of c-KIT kinase with IC₅₀s of 4 nM and 8 nM for c-Kit (wt) and c-Kit (T670I).