
STING agonist-4

产品编号: D52204

CAS: 2138300-40-8

分子式: C₃₄H₃₈N₁₂O₄

纯度: 98%

InChi Key: InChIKey=ICZSAXDKFXTSGL-UHFFFAOYSA-N

Smiles: CCn1nc(C)cc1C(Nc1nc(cc(cc2)C(N)=O)c2n1CCCCn1c(NC(c2cc(C)nn2CC)=O)nc2c1ccc(C(N)=O)c2)=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: STING

溶解性: DMSO: 10 mg/mL (14.73 mM; Need ultrasonic); H₂O: < 0.1 mg/mL (insoluble)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: STING agonist-4 is an agonist of stimulator of Interferon Genes (STING) receptor (IC₅₀: 20 nM).