
YLF-466D

产品编号: D52217

CAS: 1273323-67-3

分子式: C₂₉H₂₀ClNO₃

纯度: 98.04%

InChi Key: InChIKey=BPBOVHROVFJFAH-CYYJNZCTSA-N

Smiles: C(=C/C1C=2C(N(CC3=CC(C(O)=O)=CC=C3)C1=O)=CC=CC2)/\C4=CC=C(Cl)C=C4)/C5=CC=CC=C5

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling;Chromatin/Epigenetic

作用靶点: AMPK

溶解性: DMSO: 90mg/mL (193.2mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: YLF-466D is a newly developed AMPK activator, which inhibits platelet aggregation.