

---

## SYN1143

产品编号: D52234

CAS: 913376-84-8

分子式: C<sub>31</sub>H<sub>29</sub>FN<sub>4</sub>O<sub>5</sub>

纯度: 98.86%

InChi Key: InChIKey=UYMSIPINLJNNOU-UHFFFAOYSA-N

Smiles: CC(C)(CN(C(C)=C(C(Nc(cc1)cc(F)c1O)c1ccnc2cc(OC)ccc12)=O)C1=O)N1c1cccc1)O

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR

溶解性: DMSO:90 mg/mL (161.7 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AMG-1 is a potent inhibitor of human RON and c-Met. In vitro kinase assays showed that compound I is a potent inhibitor of human RON and c-Met with IC<sub>50</sub>s of 9 and 4 nmol/L, respectively.