
AMPA receptor modulator-2

产品编号: D52236

CAS: 2034181-36-5

分子式: C₁₅H₆F₆N₄O₅

纯度: 99.62%

InChi Key: InChiKey=BPZAQDCXSKBMKM-UHFFFAOYSA-N

Smiles: O=C1Sc(cc(cc2)-n3c4nc(C(F)(F)F)ccc4nc3C(F)(F)F)c2N1

外观: 固体粉末

作用通路: Neuroscience

作用靶点: GluR

溶解性: DMSO : 250 mg/mL (618.37 mM; Need ultrasonic); H₂O : < 0.1 mg/mL (insoluble)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AMPA receptor modulator-2 is a modulator of AMPA receptor. The pIC₅₀ value for TARPγ2 dependent AMPA receptor is 10.1.