
Glucocorticoid receptor agonist-1

产品编号: D52262

CAS: 2166375-82-0

分子式: C₃₅H₃₉NO₆

纯度: 99.51%

InChi Key: InChIKey=MPEPSOPXQDGEHP-WWBNNTNKSA-N

Smiles: C(CO)(=O)[C@]12[C@]3(C)[C@@](C[C@]1(O[C@H](O2)C4=CC=C(CC5=CC(N)=CC=C5)C=C4)[H])([C@]6([C@]([C@@H](O)C3)([C@]7(C)C(CC6)=CC(=O)C=C7)[H])[H])[H]

外观: 固体粉末

作用通路: Endocrinology/Hormones

作用靶点: Glucocorticoid Receptor

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Glucocorticoid receptor agonist-1 is a potent glucocorticoid receptor agonist with an IC₅₀ of 2.8 nM extracted from patent WO2017210471A1, compound 41[1].