
MM 11253

产品编号: D52301

CAS: 345952-44-5

分子式: C₂₈H₃₀O₂S₂

纯度: 98.62%

InChi Key: InChIKey=DPVGXUNWCVEIGT-UHFFFAOYSA-N

Smiles: O=C(O)C=1C=CC2=CC(=CC=C2C1)C3(SCCS3)C4=CC=C5C(=C4)C(C)(C)CCC5(C)C

外观: 固体粉末

作用通路: Metabolism

作用靶点: Retinoid Receptor

溶解性: DMSO : 100 mg/mL (216.1 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MM 11253 is a RAR γ antagonist with IC₅₀ of 44nM.