
BAY-85-8501

产品编号: D52307

CAS: 1161921-82-9

分子式: C₂₂H₁₇F₃N₄O₃S

纯度: 98.63%

InChi Key: InChIKey=YAJWYFPMASPAMM-HXUWFJFHSA-N

Smiles: CC(N(c1cc(C(F)(F)F)ccc1)C(N(C)[C@@H]1c(ccc(C#N)c2)c2S(C)(=O)=O)=O)=C1C#N

外观: 固体粉末

作用通路: Proteases/Proteasome

作用靶点: Serine Protease

溶解性: DMSO:200 mg/mL (421.53 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BAY-85-8501 is a selective, reversible and potent inhibitor of Human Neutrophil Elastase (HNE; IC₅₀ of 65 pM)