
GW406108X(Z/E)

产品编号: D52374

CAS: 265098-01-9

分子式: C₂₀H₁₁Cl₂NO₄

纯度: 99.71%

InChi Key: InChIKey=XKTUKRBLWOHYIL-UHFFFAOYSA-N

Smiles: O=C1NC2=C(C=C(C(C3=CC=CO3)=O)C=C2)C1=CC4=CC(Cl)=C(O)C(Cl)=C4

外观: 固体粉末

作用通路: Autophagy; Cytoskeletal Signaling

作用靶点: Autophagy; Kinesin

溶解性: DMSO: 30 mg/ml (74.96 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GW406108X(Z/E) is a mixture of different configurations of GW406108X, which is an inhibitor of Kinesin-12 and ULK1 .