
Bupranolol

产品编号: D52381

CAS: 14556-46-8

分子式: C₁₄H₂₂ClNO₂

纯度: 95.85%

InChi Key: InChIKey=HQIRNZOQPUAHHV-UHFFFAOYSA-N

Smiles: Cc1ccc(Cl)c(OCC(O)CNC(C)(C)C)c1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Adrenergic Receptor

溶解性: DMSO: 25 mg/mL (91.99 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bupranolol is an antagonist of β ₂-adrenergic receptor, and has been used for angina pectoris, hypertension, glaucoma, cardiac arrhythmia and as an antithrombotic.