
Mc-VC-PAB-SN38

产品编号: D52387

CAS: 1801838-28-7

分子式: C₅₁H₅₈N₈O₁₃

纯度: 98.74%

InChi Key: InChIKey=CNDJKXYZKSJFNL-XLPKKHSNSA-N

Smiles: O=C1N2C(C=3C(C2)=C(CC)C=4C(N3)=CC=C(O)C4)=CC5=C1COC(=O)[C@]5(OC(OCC6=CC=C(NC([C@@H](NC([C@@H](NC(CCCCN7C(=O)C=CC7=O)=O)[C@@H](C)C)=O)CCCNC(N)=O)=O)C=C6)=O)CC

外观: 固体粉末

作用通路: Antibody-drug Conjugate/ADC Related

作用靶点: Antibody-Drug Conjugates (ADCs)

溶解性: DMSO:20 mg/mL (20.18 mM,Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Mc-VC-PAB-SN38 consists of a cleavable ADC linker (Mc-VC-PAB) and SN38, which is part of an antibody-coupled reactive molecule commonly used to synthesize antibody-coupled reactive molecules (ADCs) that target sites of action.