

---

## L-Argininamide dihydrochloride

产品编号: D52406

CAS: 14975-30-5

分子式: C<sub>6</sub>H<sub>17</sub>Cl<sub>2</sub>N<sub>5</sub>O

纯度: 100%

InChi Key: InChiKey=LYMQLFYWIDCFLC-FHNDMYTFSA-N

Smiles: Cl.Cl.N[C@@H](CCCNC(N)=N)C(N)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 50 mg/mL (203.14 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: L-Argininamide (L-Arm) is used as a model compound in studies of the physicochemical characteristic of ligand binding DNA aptamers and their potential development as fluorescent aptasensors.