
Bredinin aglycone

产品编号: D52417

CAS: 56973-26-3

分子式: C₄H₅N₃O₂

纯度: 99.88%

InChi Key: InChIKey=UEWSIIBPZOBMBL-UHFFFAOYSA-N

Smiles: NC(c(nc[nH]1)c1O)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 21.67 mg/ml (170.47 mM); H₂O: 1 mg/mL (7.87 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bredinin aglycone is a purine nucleotide analog. It can be used to examine the efficiency of catalysts for the preparation of purine nucleotide analogs.