
R59949

产品编号: D52434

CAS: 120166-69-0

分子式: C₂₈H₂₅F₂N₃O₃S

纯度: 98.82%

InChi Key: InChIKey=ZCNBZFRECRPCKU-UHFFFAOYSA-N

Smiles: O=C1N(CCN2CC/C(CC2)=C(C3=CC=C(F)C=C3)\C4=CC=C(F)C=C4)C(NC5=C1C=CC=C5)=S

外观: 固体粉末

作用通路: Chromatin/Epigenetic; Cytoskeletal Signaling

作用靶点: PKC

溶解性: DMSO: 60 mg/ml (122.55 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: R 59949 is used as a Diacylglycerol kinase (DGK) inhibitor.