
LED209

产品编号: D52446

CAS: 245342-14-7

分子式: C₁₉H₁₇N₃O₂S₂

纯度: 98%

InChi Key: InChIKey=HNDRSTUKPCLQLT-UHFFFAOYSA-N

Smiles: O=S(=O)(Nc1ccccc1)c1ccc(NC(=S)Nc2ccccc2)cc1

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Antibacterial

溶解性: DMSO:10 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LED209, a specific inhibitor of quorum sensor QseC of GNB, can generate a multifunctional agent PAMAM-LED209.