
AbBy Fluor[®]680, NHS Ester

产品编号: D-9711

保存条件: Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

产品介绍: 产品参数

Ex/Em: 680/701 nm

贮存条件: -20°C避光保存

保质期: 至少可以储存6个月以上

分子量: ~1012

消光系数: 140,000